

BALOANE INDICE IOD

BALON INDICE IOD

- Sticla borosilicata 3,3 mm, termorezistenta
- Slif cu dop rodat



Cod	Produs	Mod de ambalare
BMSE82.0001	Balon indice iod, 50 ml	bucata
BMDE83.0001	Balon indice iod, 100 ml	bucata
BMNE84.0001	Balon indice iod, 200 ml	bucata
BMTE85.0001	Balon indice iod, 250 ml	bucata
BMVE86.0001	Balon indice iod, 300 ml	bucata
BMPE87.0001	Balon indice iod, 500 ml	bucata
BMOE88.0001	Balon indice iod, 1000 ml	bucata