

## PALNII SEPARARE

### PALNII SEPARARE

- Sticla borosilicata de 3,3 mm termorezistenta



### Palnie de separare cilindrica, negradata, cu robinet din sticla

Cod	Produs	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMIE316.0001	Palnie separare, cilindrica, 60 ml	9 x 80	220	bucata
PMLE317.0001	Palnie separare, cilindrica, 125 ml	9 x 80	260	bucata
PMHE77.0001	Palnie separare, cilindrica, 250 ml	10 x 110	290	bucata
PMVE78.0001	Palnie separare, cilindrica, 500 ml	10 x 110	330	bucata
PMGE79.0001	Palnie separare, cilindrica, 1000 ml	12 x 110	400	bucata

### Palnie de separare cilindrica, negradata, cu robinet din teflon

Cod	Produs	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMBE318.0001	Palnie separare, cilindrica, 60 ml	9 x 80	220	bucata
PMDE319.0001	Palnie separare, cilindrica, 125 ml	9 x 80	260	bucata
PMCE320.0001	Palnie separare, cilindrica, 250 ml	10 x 110	290	bucata
PMXE321.0001	Palnie separare, cilindrica, 500 ml	10 x 110	330	bucata
PMGE322.0001	Palnie separare, cilindrica, 1000 ml	12 x 110	400	bucata

### Palnie de separare conica, negradata, cu robinet din sticla

Cod	Produs	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMPE323.0001	Palnie separare, conica, 60 ml	10 x 70	220	bucata
PMTE324.0001	Palnie separare, conica, 125 ml	10 x 70	260	bucata
PMVE325.0001	Palnie separare, conica, 250 ml	11 x 70	290	bucata
PMTE74.0001	Palnie separare, conica, 500 ml	12 x 70	330	bucata
PMFE326.0001	Palnie separare, conica, 1000 ml	12 x 70	400	bucata
PMHE509.0001	Palnie separare, conica, 2000 ml	12 x 70	450	bucata

### Palnie de separare, conica, negradata cu robinet din teflon

Cod	Produs	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMQE327.0001	Palnie separare, conica, 60 ml	10 x 70	220	bucata
PMSE328.0001	Palnie separare, conica, 125 ml	10 x 70	260	bucata
PMXE329.0001	Palnie separare, conica, 250 ml	11 x 70	290	bucata
PMKE330.0001	Palnie separare, conica, 500 ml	12 x 70	330	bucata
PMZE331.0001	Palnie separare, conica, 1000 ml	12 x 70	400	bucata



### Palnie de separare cilindrica, gradata (gradatie orientativa), cu robinet din sticla

Cod	Produs	Gradatie (ml)	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMRE332.0001	Palnie separare, cilindrica, 60 ml	5	9 x 80	220	bucata
PMCE333.0001	Palnie separare, cilindrica, 125 ml	5	9 x 80	260	bucata

PMSE334.0001	Palnie separare, cilindrica, 250 ml	10	10 x 110	290	bucata
PMLE335.0001	Palnie separare, cilindrica, 500 ml	20	10 x 110	330	bucata
PMGE336.0001	Palnie separare, cilindrica, 1000 ml	50	12 x 110	400	bucata

**Palnie de separare cilindrica, gradata (gradatie orientativa), cu robinet din teflon**

Cod	Produce	Gradatie (ml)	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMTE337.0001	Palnie separare, cilindrica, 60 ml	5	9 x 80	220	bucata
PMHE338.0001	Palnie separare, cilindrica, 125 ml	5	9 x 80	260	bucata
PMJE339.0001	Palnie separare, cilindrica, 250 ml	10	10 x 110	290	bucata
PMWE340.0001	Palnie separare, cilindrica, 500 ml	20	10 x 110	330	bucata
PMYE341.0001	Palnie separare, cilindrica, 1000 ml	50	12 x 110	400	bucata

**Palnie de separare conica, gradata (gradatie orientativa), cu robinet din sticla**

Cod	Produce	Gradatie (ml)	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMPE342.0001	Palnie separare, conica, 60 ml	5	10 x 70	220	bucata
PMNE343.0001	Palnie separare, conica, 125 ml	5	10 x 70	260	bucata
PMXE344.0001	Palnie separare, conica, 250 ml	10	11 x 70	290	bucata
PMME345.0001	Palnie separare, conica, 500 ml	20	12 x 70	330	bucata
PMVE346.0001	Palnie separare, conica, 1000 ml	50	12 x 70	400	bucata
PMLE508.0001	Palnie separare, conica, 2000 ml	50	12 x 70	450	bucata

**Palnie de separare conica, gradata (gradatie orientativa), cu robinet din teflon**

Cod	Produce	Gradatie (ml)	Ø x L tija [orientativ] (mm)	Lungime totala [orientativ] (mm)	Mod de ambalare
PMLE347.0001	Palnie separare, conica, 60 ml	5	10 x 70	220	bucata
PMGE348.0001	Palnie separare, conica, 125 ml	5	10 x 70	260	bucata

PMBE349.0001	Palnie separare, conica, 250 ml	10	11 x 70	290	bucata
PMZE350.0001	Palnie separare, conica, 500 ml	20	12 x 70	330	bucata
PMRE351.0001	Palnie separare, conica, 1000 ml	50	12 x 70	400	bucata